

Bacterial Ident 5 1 BDC 0 Tc 0 Tc ca 0107 Tn

Performance

Days Performed

Daily; Monday – Sunday

Report Available (TAT) – (Once received at MDL)

2 – 4 days

Specimen Retention Time

7 days

Method Description

Identification methods may include any of the following: conventional biochemical testing, matrix-assisted laser desorption/ionization time-of-flight (MALDI-TOF) mass spectrometry, and commercial identification panels. Genus and species are reported on aerobic bacterial isolates whenever possible.